

Flow cytometry solutions

## The MACSQuant® VYB



Excite the color that counts

# Achieve pure analytical power

The yellow laser expands your possibilities

## The MACSQuant® VYB is designed for maximal versatility

The MACSQuant® VYB delivers all the performance, automation, compactness, and convenience of the MACSQuant Analyzer but with a uniquely configured optical layout, featuring violet, yellow, and blue lasers. With this new optical layout in combination with ten optical detection channels, the MACSQuant VYB is a powerful and versatile flow cytometer for virtually every lab.

### A richer dataset in less time

By incorporating a yellow 561 nm laser and ten optical detection channels, the MACSQuant VYB offers researchers a greater choice of fluorescent proteins and conjugates to perform sophisticated multicolor experiments.

- Analyze multiple fluorescent proteins at once, such as mCherry, DsRed, GFP, CFP, and many more.
- Profit from optimal PE sensitivity for enhanced detection of even PE<sup>dim</sup> populations.
- Rely on an optical configuration optimized for the detection of all commonly used fluorochromes.
- Reliably detect APC signals and APC tandem dyes, such as APC-Vio770™, APC-Cy™7, or APC-H7 with the powerful 100 mW yellow laser.
- Obtain comparable stain indices and mean fluorescence intensities compared with red 635 nm laser systems.

Excitation laser	Channel (filter)	Antibody/parameter	Fluorescent protein
Violet 405 nm	V1 (450/50 nm)	VioBlue®	CFP
	V2 (525/50 nm)	VioGreen™	vGFP, GFP
Blue 488 nm	B1 (525/50 nm)	FITC	GFP, YFP
	B2 (614/50 nm)	PI	Lss-mKate
Yellow 561 nm	Y1 (586/15 nm)	PE	
	Y2 (615/20 nm)	PE-Texas Red®	mCherry, DsRed, tdTomato
	Y3 (661/20 nm)	PE-Cy5/APC	mKate
	Y4 (750 nm/LP)	PE-Vio770/APC-Vio770	mPlum
Yellow 561 nm	FSC (561/10 nm)	Forward scatter	
	SSC (561/10 nm)	Side scatter	

**Table 1:** Optical configuration of the MACSQuant VYB.



**"These results [...] demonstrate the ability of the MACSQuant VYB to resolve multiple fluorescent proteins and organic dyes at the same time."**


Kueh, H. Y. and Diamond, R. (2012) MACS&more 14–1: 8–10.

### Minimal compensation

Combine fluorescent protein reporters with antibody surface stainings with ease. With PE detection from the yellow 561 nm laser and FITC or GFP detection from the blue 488 nm laser, there is little need for compensation with the MACSQuant VYB.

### Optimal PE performance

PE detection from the yellow 561 nm laser means greater sensitivity for the discrimination of dim populations and improved detection with PE tandems such as PE-Vio770 or PE-Cy5.



**Figure 1:** Cells expressing GFP under an Oct4 promoter and stained for Anti-Feeder-PE, a dim marker for fibroblasts, were analyzed using the MACSQuant VYB. Excitation of both PE and GFP with the blue 488 nm laser only minimally resolves the two signals (**A**). In comparison, excitation of PE by the yellow 561 nm laser and of GFP by the blue 488 nm laser fully resolves these two signals without the need for any compensations (**B**).


# Discover greater versatility

Multicolor analysis broadens your research horizons

## Achieve new answers

### Discriminate APC populations with ease


With four-laser functionality in a three-laser instrument, there are no constraints in your choice of fluorescent proteins and conjugates.



**Figure 2:** Lymphocytes were stained with CD3-FITC, CD4-PE, and CD8-APC. The MACSQuant Analyzer 10 uses the blue 488 nm laser for excitation of FITC and PE, and the red 635 nm laser for APC (A–C), while the MACSQuant VYB uses the blue 488 nm laser for FITC only, and the yellow 561 nm laser for PE and APC (D–F). APC resolution with the yellow 561 nm laser is equal to the resolution from a red 635 nm laser. MFI: mean fluorescence intensity; SI: stain index.

### Classical immunophenotyping

The MACSQuant® VYB allows researchers to perform elaborate multiplex investigations by detecting a greater number of colors simultaneously. This leads to the acquisition of more meaningful data in less time.



**Figure 3:** Classic immunophenotyping of peripheral blood stained with a panel of six standard dyes and antibodies (B–F): propidium iodide (PI), CD45-VioBlue, CD3-FITC, CD4-PE, CD8-APC, and CD45RA-APC-Cy7.

## Rare cell analysis without the waiting

Save time without compromising on the quality and sensitivity of your rare cell analyses. With the MACSQuant Instruments you can easily enrich rare cell populations prior to multiparametric cell analysis.

- Significantly reduce lengthy acquisition times.
- Deliver highly sensitive rare cell analysis.
- Allow easy analysis of particularly challenging cell samples, such as rare cell populations in high volumes of buffer, or precious low volume samples.

Product	Order no.
<b>Instruments</b>	
MACSQuant VYB	130-096-116
MACSQuant Analyzer	130-092-197
MACSQuant Analyzer 10	130-096-343
MACS MiniSampler	130-092-535
<b>Service contracts</b>	
MACSQuant VYB Full Service L	160-001-637
MACSQuant Full Service	160-001-286
<b>Accessories</b>	
Chill 5 Rack	130-092-951
Chill 96 Rack	130-094-459
MACSQuant Columns	130-094-458
MACSQuant Calibration Beads	130-093-607
MACSQuant Running Buffer	130-092-747
MACSQuant Washing Solution	130-092-749
MACSQuant Storage Solution	130-092-748
MACSQuant Washing & Storage Solution Kit	130-092-801
MACS Bleach Solution	130-093-663
MACSQuant Starting Buffer Kit	130-094-190

► [macsquant.com](http://macsquant.com)



Miltenyi Biotec

**Germany/Austria/  
Switzerland**  
Miltenyi Biotec GmbH  
Friedrich-Ebert-Straße 68  
51429 Bergisch Gladbach  
Germany  
Phone +49 2204 8306-0  
Fax +49 2204 85197  
macs@miltenyibiotec.de

**USA/Canada**  
Miltenyi Biotec Inc.  
2303 Lindbergh Street  
Auburn, CA 95602, USA  
Phone 800 FOR MACS  
Phone +1 530 888 8871  
Fax +1 530 888 8925  
macs@miltenyibiotec.com

**Australia**  
Miltenyi Biotec  
Australia Pty. Ltd.  
Unit 16A, 2 Eden Park Drive  
Macquarie Park NSW 2113  
Australia  
Phone +61 2 8877 7400  
Fax +61 2 9889 5044  
macs@miltenyibiotec.com.au

**Benelux**  
Miltenyi Biotec B.V.  
Schipholweg 68 H  
2316 XE Leiden  
The Netherlands  
macs@miltenyibiotec.nl  
**Customer service**  
**The Netherlands**  
Phone 0800 4020120  
Fax 0800 4020100  
**Customer service Belgium**  
Phone 0800 94016  
Fax 0800 99626  
**Customer service Luxembourg**  
Phone 800 24971  
Fax 800 24984

**China**  
Miltenyi Biotec GmbH  
Shanghai Office  
Rm. 2309  
No. 319 Xianxia Rd.  
Shanghai 200051, P.R. China  
Phone +86 21 62351005  
Fax +86 21 62350953  
macs@miltenyibiotec.com.cn

**France**  
Miltenyi Biotec SAS  
10 rue Mercoeur  
75011 Paris, France  
Phone +33 1 56 98 16 16  
Fax +33 1 56 98 16 17  
macs@miltenyibiotec.fr  
**Italy**  
Miltenyi Biotec S.r.l.  
Via Persicetana, 2/D  
40012 Calderara di Reno (BO)  
Italy  
Phone +39 051 6 460 411  
Fax +39 051 6 460 499  
macs@miltenyibiotec.it

**Japan**  
Miltenyi Biotec K.K.  
Nittsu-Eitai Building 5F  
16-10 Fuyuki, Koto-ku,  
Tokyo 135-0041, Japan  
Phone +81 3 5646 8910  
Fax +81 3 5646 8911  
macs@miltenyibiotec.jp

**Nordics and Baltics**  
Miltenyi Biotec Norden AB  
Scheelevägen 17  
223 70 Lund  
Sweden  
macs@miltenyibiotec.se  
**Customer service Sweden**  
Phone 0200-111 800  
Fax 046-280 72 99  
**Customer service Denmark**  
Phone 80 20 30 10  
Fax +46 46 280 72 99  
**Customer service**  
**Norway, Finland, Iceland,  
and Baltic countries**  
Phone +46 46 280 72 80  
Fax +46 46 280 72 99

**Singapore**  
Miltenyi Biotec Asia Pacific Pte Ltd.  
100 Beach Road  
#28-06 to 28-08 Shaw Tower  
Singapore 189702  
Phone +65 6238 8183  
Fax +65 6238 0302  
macs@miltenyibiotec.com.sg

**South Korea**  
Miltenyi Biotec Korea  
Petra building  
659 Yeoksam-dong  
Gangnam-gu  
Seoul, South Korea  
Phone +82 2 555 1988  
Fax +82 2 555 8890  
macs@miltenyibiotec.co.kr

**Spain**  
Miltenyi Biotec S.L.  
C/Luis Buñuel 2  
Ciudad de la Imagen  
28223 Pozuelo de Alarcón (Madrid)  
Spain  
Phone +34 91 512 12 90  
Fax +34 91 512 12 91  
macs@miltenyibiotec.es  
**United Kingdom**  
Miltenyi Biotec Ltd.  
Almac House, Church Lane  
Bisley, Surrey GU24 9DR, UK  
Phone +44 1483 799 800  
Fax +44 1483 799 811  
macs@miltenyibiotec.co.uk

[www.miltenyibiotec.com](http://www.miltenyibiotec.com)

Miltenyi Biotec provides products and services worldwide. Visit [www.miltenyibiotec.com/local](http://www.miltenyibiotec.com/local) to find your nearest Miltenyi Biotec contact.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS, MACSQuant, Vio770, VioBlue, and VioGreen are either trademarks or registered trademarks of Miltenyi Biotec GmbH. All other trademarks mentioned in this document are the property of their respective owners and are used for identification purposes only. Copyright © 2013 Miltenyi Biotec GmbH. All rights reserved.