

1. Description

1.1 Background information

1.2 Technical specifications

1.3 Further information

1.4 Applications

1. Description

This product is for research use only.

Components 6-well or 24-well or 96-well plates, black, silicon foil bottom, with lid, sterile and single packed, 4 pieces per box.

Product	Order no.
Gas-Permeable Culture Plate (6 well)	150-000-363
Gas-Permeable Culture Plate (24 well)	150-000-362
Gas-Permeable Culture Plate (96 well)	150-000-364

Storage Store Gas-Permeable Culture Plates dry and protected from light at room temperature (19–25 °C).
The expiration date is indicated on the box label.
Do not use after this date.

1.1 Background information

Gas-Permeable Culture Plates have been designed for cell-based assays with special requirements for control of gas partial pressure in the cellular environment. Applications are high density cell cultures, cultivation of cells with high oxygen demands or in exactly controlled oxygen environments like hypoxic experimental setting. The plate bottom is made from thin silicone foil (170 µm) and enables rapid gas exchange between cell culture media and atmosphere.

1.2 Technical specifications

- Plate material: polystyrene, black
- Bottom material: silicone foil, transparent
- Lid material: polystyrene, transparent
- Recommended working temperature between 4 °C and 40 °C
- Recommended working volume per well: 3–4 mL (6 well), 500–1000 µL (24 well), 100–200 µL (96 well)
- Dimensions (ANSI/SLAS guideline compliant):

	6-well plate	24-well plate	96-well plate
Length [mm]	127.76	127.76	127.76
Width [mm]	85.48	85.48	85.48
Height [mm]	15	15	14.35
Well diameter [mm]	32	13	6
Well-well distance x/y [mm]	36	18	9
Distance well center A1 to short side [mm]	27.88	18.88	14.38
Distance well center A1 to long side [mm]	24.74	15.74	11.24
Distance bottom to focal plane [mm]	0.4	0.4	0.4
Bottom thickness [mm]	0.17	0.17	0.17

1.3 Further information

- Gas-Permeable Culture Plates ensure well-controlled humidity >95% in the incubator to prevent accelerated evaporation. Spacer plates might be used to ensure a gap between the plate bottom and the incubator shelf.
- Best performance is ensured when the plates are not stacked (ensure bottom ventilation).
- Recommended for suspension cells. Adherence-dependent cells need biological coatings appropriate for the used cell type.

1.4 Applications

Gas-Permeable Culture Plates have been designed for cell-based assays, such as:

- T cell activation
- NK cell activation
- Controlled hypoxia experiments

Refer to www.miltenyibiotec.com for all data sheets and protocols. Miltenyi Biotec provides technical support worldwide. Visit www.miltenyibiotec.com/local to find your nearest Miltenyi Biotec contact.

Legal notices

Limited product warranty

Miltenyi Biotec B.V. & Co. KG and/or its affiliate(s) warrant this product to be free from material defects in workmanship and materials and to conform substantially with Miltenyi Biotec's published specifications for the product at the time of order, under normal use and conditions in accordance with its applicable documentation, for a period beginning on the date of delivery of the product by Miltenyi Biotec or its authorized distributor and ending on the expiration date of the product's applicable shelf life stated on the product label, packaging or documentation (as applicable) or, in the absence thereof, ONE (1) YEAR from date of delivery ("Product Warranty"). Miltenyi Biotec's Product Warranty is provided subject to the warranty terms as set forth in Miltenyi Biotec's General Terms and Conditions for the Sale of Products and Services available on Miltenyi Biotec's website at www.miltenyibiotec.com, as in effect at the time of order ("Product Warranty"). Additional terms may apply. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING IF A PRODUCT IS SUITABLE FOR CUSTOMER'S PARTICULAR PURPOSE AND APPLICATION METHODS.

Technical information

The technical information, data, protocols, and other statements provided by Miltenyi Biotec in this document are based on information, tests, or experience which Miltenyi Biotec believes to be reliable, but the accuracy or completeness of such information is not guaranteed. Such technical information and data are intended for persons with knowledge and technical skills sufficient to assess and apply their own informed judgment to the information. Miltenyi Biotec shall not be liable for any technical or editorial errors or omissions contained herein.

All information and specifications are subject to change without prior notice. Please contact Miltenyi Biotec Technical Support or visit www.miltenyibiotec.com for the most up-to-date information on Miltenyi Biotec products.

Licenses

This product and/or its use may be covered by one or more pending or issued patents and/or may have certain limitations. Certain uses may be excluded by separate terms and conditions. Please contact your local Miltenyi Biotec representative or visit Miltenyi Biotec's website at www.miltenyibiotec.com for more information.

The purchase of this product conveys to the customer the non-transferable right to use the purchased amount of the product in research conducted by the customer (whether the customer is an academic or for-profit entity). This product may not be further sold. Additional terms and conditions (including the terms of a Limited Use Label License) may apply.

CUSTOMER'S USE OF THIS PRODUCT MAY REQUIRE ADDITIONAL LICENSES DEPENDING ON THE SPECIFIC APPLICATION. THE CUSTOMER IS SOLELY RESPONSIBLE FOR DETERMINING FOR ITSELF WHETHER IT HAS ALL APPROPRIATE LICENSES IN PLACE. Miltenyi Biotec provides no warranty that customer's use of this product does not and will not infringe intellectual property rights owned by a third party. BY USE OF THIS PRODUCT, THE CUSTOMER AGREES TO BE BOUND BY THESE TERMS.

Trademarks

The Miltenyi Biotec logo is a registered trademark or trademark of Miltenyi Biotec and/or its affiliates in various countries worldwide. All other trademarks mentioned in this publication are the property of their respective owners and are used for identification purposes only.

Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.