

Spatial information on receptor clustering in flow cytometry using FRET

Kerstin von Kolontaj¹, Gabor Horváth², Eicke Latz², and Martin Büscher¹ ¹Miltenyi Biotec B.V. & Co. KG, Bergisch Gladbach, Germany | ²Institute of Innate Immunity, University of Bonn, Bonn, Germany

Introduction

Flow cytometry is commonly used to detect the overall fluorescent signal of a whole cell and thus to analyze whether or not certain target molecules are present on the cell. However, no information on the distribution or clustering of the molecules of interest can be gained from these signals. Yet, the clustering of receptors or the interaction of proteins plays a key role in almost all biological functions.

like two-hybrid systems, are usually very laborious and highly time consuming. Foerster resonance energy transfer (FRET)¹ is often used to identify protein-protein interactions by confocal microscopy².

However, this technique requires expert knowledge and produces vast amounts of data. In order to overcome these limitations, we developed a FRET-based flow cytometry analysis

Program application: Analysis of CD3-CD4 FRET on T cells

The FRET program was applied to measure the dynamic clustering of CD3 and CD4 receptors on the cell surface after T cell activation. Here, PBMC were isolated from buffy coat and the T cells were activated using 1 µg staphylococcal enterotoxin B (SEB). (A) The analysis using the FRET program showed that the FRET efficiency significantly increased after a 2-minute incubation with SEB and decreased again 5 hours after activation.

(B) The clustering could be visualized using

equally on the cell surface. After a 2-minute activation with SEB, small activation clusters of CD3 and CD4, the so-called microclusters, began to form. After five hours of incubation, large clusters of CD3 had formed⁴. (C) The interaction between CD3 and CD4 was validated using a colocalization analysis tool for CM. The normalized colocalization data generated by CM are highly similar to the data resulting from the FRET analysis program. However, it is obvious that for non-activated cells the colocalization coefficient is significantly higher than the FRET efficiency, which is due to the lower resolution of the confocal microscope.

Therefore, plenty of research efforts are spent in this field. Conventional investigation methods,

program, which automatically provides spatial information about protein interactions on the $cell^3$.

Methods

Automatic FRET measurement

The FRET Express Mode Program automatically measures and determines the FRET efficiency on a cell-by-cell basis on the MACSQuant[®] Analyzer 10.

Four cell samples are required to accurately determine FRET efficiency:

- Blank 1. Subtraction of background fluorescence
- 2./3. Donor only and Acceptor only Determination of cross-talk between the donor and acceptor channels
- FRET sample 4. Calculation of FRET efficiency in the sample

The controls and FRET samples are measured and gated for the target populations by the program. Briefly, for the determination of the FRET efficiency, relative changes in donor and acceptor fluorescence intensities are registered. Based on this and taking into account the results obtained from the control samples, the rate of transferred energy is then calculated accurately by the FRET Express Mode Program in a fully automated manner.

confocal microscopy (CM). Cells were stained with CD3-FITC and CD4-APC. For the non-activated cells, CD3 and CD4 receptors were distributed Α *p < 0.01 **p < 0.01









No activation

2 min activation

5 h activation

Figure 3

В





Figure 1

Results

Validation of program functionality

We validated the sensitivity and accuracy of the FRET program using FRET standard beads with defined linker lengths between the fluorochromes. FRET efficiency readouts generated automatically by the program were compared to established non-automated calculations. The high correlation FRET between the two methods (R2 = 0.95) proved the functionality of the FRET program.



The FRET program enabled us to analyze homoclustering of the CD3 receptor after T cell activation. A single epitope of CD3 was stained with either CD3-VioBlue[®] or CD3-FITC, and cells were activated with SEB for 2 minutes. T cell activation led to a significant increase in FRET

efficiency due to CD3 cluster formation⁴. T cell activation could be inhibited by pre-incubating SEB with an anti-SEB blocking antibody (A) before adding SEB to the cells and was reduced by adding SEB and its blocking antibody simultaneously (B).





Figure 4

Α

Conclusion

We developed a program which enables the automated flow cytometric measurement of FRET in living cells. This so-called FRET Express Mode program allows for the detection of proteinprotein interactions on cells. The program

References

- 1. Förster, T. (1948) Annals of Physics 2: 55–75.
- 2. Nagy, P. et al. (1998) Eur. Biophys. 27: 377–389.
- 3. Stryer, L. et al. (1978) Annu. Rev. Biochem. 47:819-846.

Unless otherwise specifically indicated, Miltenyi Biotec products and services are for research use only and not for therapeutic or diagnostic use. MACS, MACSQuant, and VioBlue are registered trademarks of Miltenyi Biotec and/or its affiliates in various countries worldwide. Copyright © 2020 Miltenyi Biotec and/or its affiliates. All rights reserved.

overcomes certain limitations of conventional flow cytometry, as it not only provides information on whether or not a particular molecule is present on a cell, but also automatically analyzes its spatial distribution and interaction with other molecules.

Therefore, the FRET Express Mode Program will simplify the analysis of biologically relevant changes within a cell.

4. Huppa, J.B. et al. (2010) Nature 18: 963–967.

universität**bonn**